

Sterile adhesive plaster for covering the site of puncture

ELEBAN PRESTAT

MEDICAL DEVICE



ELEBAN PRESTAT is a sterile adhesive plaster with a pad, designed for covering the site of puncture. It can be used for blood stanching. (homostatic material)



● Adhesive tape

Excellent at elasticity and the strength of oppression is adjustable

Urethane-made nonwoven sheet with low-allergic acrylic adhesive

● Wrapping

Packed individually, sterilized with EO

● Absorbent pad

The pad made of cotton expands with absorbing blood. That increases the effect of blood stanching



size: 35mm × 80mm
(pad: 20mm in diameter
3.7mm in thickness)

Cotton-made compressed pad with right hardness

Hakuzo Medical Corporation

4-9, 2 chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan
Tel; 81-6-6809-4319 Fax; 81-6-6945-0539
E-mail; trade@hakuzo.co.jp

MADE IN JAPAN